

Accessories

SURGICAL LASER ACCESSORIES

Superior accessories for successful procedures



COMPETENCE IN SURGICAL LASER SINCE 1989
LASER TECHNOLOGY - MADE IN GERMANY
WWW.LISALASER.DE



LISA laser products provides various accessories and helpful devices for many laser applications in different medical specialities and complementary product lines.

For our range of fibres please see the Surgical Laser Fibres brochure.

For uretherorenoscopes, nephroscopes and complementary instruments please see the ScopeLine brochures.

Fibre Insertion and Fibre Fixation



TuohyBorst

Port Seal for fixation of laser fibres to endoscopes and laser instruments, connects to female Luer Lock, sterile on delivery, single use

1 package á 25 pieces

200 001 030



FlexGuard

laser fibre insertion sheath, sterile on delivery, for single use, outer diameter 0.9 mm for 3.6 F working channels, clearance 0.6 mm matching with FlexiFib/ FlexiFib-SU and LithoFib, length 1000 mm, DE 10 2005 017 204, U.S. Patent No. 7,775,969.

101 503 263

Spine Surgery



SpineLas 375-7.5

side fire applicator, for laser foraminoplasty and discectomy, working length 375 mm, outer diameter 7.5 F, includes non-detachable laser fibre, sterile on delivery, for single use

101 503 251

Eyewear

**Safety Eyewear**

Laser protective eyewear (plastic goggle fits over glasses) for Sphinx and RevOLix product families

101 503 141

**Safety Eyewear**

Colour neutral laser safety glasses for Sphinx and RevOLix product families

101 503 399

**Safety Eyewear**

Colour neutral laser safety glasses (for spectacle wearers) for Sphinx and RevOLix product families

101 503 400

Tools

**Power Meter PED**

Laser power meter

101 503 140

**WarnLED**

laser warning sign, hangs on door handles etc., with flashing LEDs

101 503 223

Warning Sign

laser warning sign for marking entrance doors

101 503 331

**FibreScope 200x**

for fibre connector inspection, magnification 200 times, illumination (not shown)

101 503 417

**FibreScope 50x**

for fibre connector inspection, magnification 50 times, w/o illumination

101 503 415

Laser Trolley
for Revolix jr. including drawer

101 510 001



IMPORTANT NOTICE:

The information provided is a general overview of potential clinical applications of the described products.

National health care regulations vary between countries and may exclude certain clinical applications at your location.

The user assumes responsibility to be updated about national deviations from the applications listed above.

In the USA the products are not intended for use in clinical applications in neurosurgery.

U.S. federal law restricts these devices to sale by or on the order of a physician

Specifications are subject to change without notice.

Made in Germany 2017-07
Brochure Part No. 036 008 213



LISA laser products OHG
Albert - Einstein - Str. 1-9
37191 Katlenburg-Lindau
Germany
fon: +49 5556 - 9938 - 0
fax: +49 5556 - 9938 - 10
info@lislaser.de www.lislaser.de



LISA LASER PRODUCTS
COMPETENCE IN SURGICAL LASER

